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
EXACTLY THE GENERATOR THAT WAS NEEDED

The design of this equipment makes it possible to supply power, using a single generator, to one or two completely separate diagnostic rooms, without having to install the control unit in the middle of the radiology rooms. The small control consoles placed in each room permit the operator to program the automatic anatomical techniques separately, according to specific operating needs. As is readily apparent, this exclusive approach not only results in significant economic saving, but also allows for the best possible use to be made of the radiology rooms


TECHNICAL SPECIFICATIONS OF THE H.F. GENERATOR

CLINIX-50 KV	CLINIX-80 KV
600 mA at 80 KV	1000 mA at 80 KV
500 mA at 100 KV	800 mA at 100 KV
400 mA at 125 KV	640 mA at 125 KV
320 mA at 150 KV	320 mA at 150 KV

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
CLINIX - HF 650

100 KHz HIGH FREQUENCY X-RAY GENERATORS
SERIES 50-80 KW ANATOMICAL PROGRAMMABLE

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PRESENTATION

The Clinix high-frequency 50-80 KV microprocessor-controlled generators is designed for all applications in the field of major radiological diagnostics. Its construction permits easy, secure maneuvering of all the controls, giving the operator the option of entering the anatomical radiological programs in the memory or using the free technique.

The extremely limited dimensions of the consoles for 1-2 tubes, designed to meet the demands of radiological diagnostics, guarantee the utmost reliability, in keeping with all the latest technological advances.

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
PRINCIPAL FEATURES

- Microprocessor control panel
- Automatic anatomical techniques
- Technique using mAs-KV
- Free technique: filament - KV - mA - sec.
- Technique with AEC exposure meter
- mAs following exposure of rays
- Pre-set mAs levels
- KV radiography from 40 to 150 in 111 steps at increments of ± 1 KV
- mA radiography from 10 to 600-1000 in 43 steps
- 1-2 work stations
- 2-4 accessories
- Control for disengagement of Potter Bucky
- Selection: 3 different body shapes
- AEC automatic exposure meter (option) - controls for: selection of field, correction of blackening, correction of sensitivity of screens
- High-speed anode rotation (option)
- KV radioscopy from 40 to 120, continuous adjustment
- mA radioscopy from 0 to 4 in 8 steps
- Radioscopy control: automatic-manual
- Display radioscopy time
- Alphabetical/numerical display for visualization of messages to operator, making it possible to monitor all the functions
- All the controls come with highly sensitive push-buttons visualized with LEDs

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- 1 Button controlling radiography timing
- 2 Button controlling radiography kv
- 3 Radiography time display
- 4 Radiography KV display
- 5 Luminous alphabetical/numerical display
- 6 mA radiography display
- 7 Button controlling radiography mA
- 8 Luminous display for anatomical regions selected
- 9 Keyboard for selection of: anatomical regions, organs and projections
- 10 Button to return alphabetical/numerical display to anatomical-region format
- 11 Button for selecting focal length
- 12 Button for selecting 1° tube - 1° potter
- 13 Button for selecting 1° tube - 2° potter
- 14 Button for selecting 2° tube - 1° potter
- 15 Button for selecting 2° tube - 2° potter
- 16 Button for engaging-disengaging potter bucky selected
- 17 Button for settings for thin, normal, stocky individuals
- 18 AEC blackening display
- 19 Button for selecting AEC blackening
- 20 Luminous signal for passage of radiography rays
- 21 Button for ON-OFF and re-set of AEC alarm
- 22 Buttons for selection of AEC screens
- 23 Button for choice of principal AEC setting
- 24 Luminous button for preparation prior to radiography ray phase
- 25 Button for radiography rays
- 26 ON-OFF buttons
- 27 Button for re-set of radioscopy timing
- 28 Button for selecting: automatic-manual radioscopy
- 29 Display of radioscopy time
- 30 Luminous signal for passage of radioscopy rays
- 31 Button for regulating radioscopy KV
- 32 Button for regulating radioscopy mA
- 33 Radioscopy KV display
- 34 Radioscopy mA display

PRINCIPAL MESSAGES TRANSCRIBED ON THE ALPHABETICAL/NUMERICAL DISPLAY:
- anatomical regions - organs - projections - thickness of patient - focal length - tube - potter - thermal units - mAs - RX ready - error-free technique.

The operator can modify the program by using a secret electronic combination.